

PCR 96 WELL PLATES COMPATIBILITY CHART

	MB-P96	*MB-P96-T	MB-Q96-LN	MB-PSM	MB-BQSM
ABI					
ProFlex 96-Wells Block	●	●		●	
ProFlex 32-Wells Block		●		●	
MiniAmp & MiniAmp Plus	●	●		●	
SimpliAmp	●	●		●	
VeritiPro 96-Wells Block	●	●		●	
Veriti 96-Wells Block (0.2mL model)	●	●		●	
Veriti Fast 96-Wells Block (0.1mL model)			●	●	
2720	●	●		●	
9700	●	●		●	
9800 Fast 96-Wells Block (0.1mL model)			●	●	
Bio-Rad					
CFX 96 & CFX 96 Connect			●		●
CFX Opus 96			●		●
T100	●	●		●	
C1000 (Touch) & S1000 (96-Wells Reaction Module)	●	●	●	●	
C1000 (Touch) & S1000 (Dual-48 Wells Reaction Module)		●		●	
iCycler (96-Wells Reaction Module)	●	●		●	
iCycler (Dual-48 Wells Reaction Module)		●		●	
MyCycler	●	●		●	
DNA Engine Tetrad 2 (96v Alpha Unit)	●	●	●	●	
DNA Engine Tetrad 2 (48-48 Alpha Unit)		●		●	

*When tearing P96-T to 8-strip tubes and using it on the ABI machine, please use the correct retainer/tray, load at least two 8-strip tubes per time and put them in adjacent columns of A3-A10 and H3-H10 to prevent deformation. For more information, please refer to: qPCR 8-Strip and 4-Strip Tubes Compatibility Chart, <https://ppt.cc/f7km5x>.