

技術資料表 Technical Data Sheet



GUNSTER BIOTECH

~Best Choice for Life Sciences~

貨號 Product NO.	MB-OPENER	發布日 Date	2024/3/6 6-March-2024
品名 Description	八連排多功能開關器 8-Strip Capping Gadget	版本 Version	第三版 Version III

產品說明 Specification				產品圖 Drawing
產品尺寸 Dimension	寬度 / 高度 Width 78.51 ±0.1mm / Height 106.46 ±0.1mm			
顏色 Color	透明 Clear			
包裝 Packaging	1 個/袋, 15 袋/盒, 20 盒/箱 1 pcs/ bag, 15 bags/ box, 20 boxes/ carton			
產品容量 Volume		儲存條件 Storage	室溫 Room temperature	
原料 Material	<ul style="list-style-type: none">● 本產品原料符合美國 FDA、聚苯乙烯● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)● 不添加殺菌劑(二 (2-羥乙基), 二甲基十二烷基銨)● 不添加增塑劑(鄰苯二甲酸鹽)● MB-OPENER is made of FDA and USP VI standard virgin Polystyrene (PS)● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMDA))● NOT INCLUDING: Plasticizers (phthalates)			
產品特色 Features	<ul style="list-style-type: none">● 適用操作環境溫度-30~70℃● 多功能設計, 可用於 8 連排管開蓋與關蓋使用, 亦可當 96 孔封膜刮板使用● 8 連排開關器之刮板處,光滑平整, 可避免對 96 孔封膜造成刮傷● 透明設計, 益於開關蓋操作● 不可高溫高壓滅菌, 可用酒精擦拭殺菌● 可用於 1.5ml 離心管● Saving more PCR operating time by using this 6-way tool● Capping and decapping separate cap on 0.1 and 0.2ml 8-Strip● Decapping attached cap on 0.2ml 8-Strip● Decapping 0.2ml PCR single tube● Decapping 1.5/1.7ml microcentrifuge tube● Scrapping adhesive film to ensure tight adhesion between film and PCR plate			

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段 192 號 9 樓之 1
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332
傳真(Fax): +886-2-2677-2339
官網(Site): www.gunster.com.tw
電郵(Email): gunster@gunster.com.tw

Eiffie Wu

Product Manager

吳佩儒
產品經理