

# 技術資料表 Technical Data Sheet



**GUNSTER BIOTECH**

~Best Choice for Life Sciences~

貨號 Product NO.	MB-E17	發布日 Date	2023/6/20 20-Jun-2023
品名 Description	1.5 ml 離心管 1.7(1.5) ml Microcentrifuge Tubes	版本 Version	第二版 Version II

產品說明 Specification				產品圖 Drawing
產品尺寸 Dimension	管口直徑 / 管身高度 Ø 9.34 ±0.01mm / Height 40.48 ±0.1mm			
顏色 Color	透明 Clear			
包裝 Packaging	500 個/袋, 10 袋/箱 500 pcs/bag, 10 bags/ transport carton			
產品容量 Volume	1.5 ml / tube	儲存條件 Storage	室溫 Room temperature	
原料 Material	<ul style="list-style-type: none"><li>● 本產品原料符合美國 FDA、美國藥典第六級認證之聚丙烯</li><li>● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)</li><li>● 不添加殺菌劑(二 (2-羥乙基)，二甲基十二烷基銨)</li><li>● 不添加增塑劑(鄰苯二甲酸鹽)</li><li>● MB-E17 is made of FDA and USP VI standard virgin Polypropylene</li><li>● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)</li><li>● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMDA))</li><li>● NOT INCLUDING: Plasticizers (phthalates)</li></ul>			
產品特色 Features	<ul style="list-style-type: none"><li>● 大面積離心蓋設計</li><li>● 最低刻度為 100 µl</li><li>● 可承受 20,000 RCF 離心力</li><li>● 管蓋和管側面具霧面設計，易標記</li><li>● 無 DNA 分解酶，RNA 分解酶，人類 DNA，PCR 抑制劑，ATP，內毒素</li><li>● 可高溫高壓滅菌(121°C, 15 psi, 15 分)</li><li>● Large-cap</li><li>● Labeled graduations at 100 µl</li><li>● Spin: 20,000 RCF</li><li>● Easy labeling on Frosted side and cap</li><li>● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin.</li><li>● Autoclavable (121°C, 15 psi, 15 minutes)</li></ul>			

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段 192 號 9 樓之 1  
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,  
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332  
傳真(Fax): +886-2-2677-2339  
官網(Site): [www.gunster.com.tw](http://www.gunster.com.tw)  
電郵(Email): [gunster@gunster.com.tw](mailto:gunster@gunster.com.tw)

*Eiffie Wu*

Product Manager

吳佩儒  
產品經理