

# 技術資料表 Technical Data Sheet



**GUNSTER BIOTECH**

~Best Choice for Life Sciences~

|                   |  |               |                          |
|-------------------|--|---------------|--------------------------|
| 貨號<br>Product NO. | MB-P02-A   | 發布日<br>Date   | 2023/6/20<br>20-Jun-2023 |
| 品名<br>Description | 0.2ml PCR 平蓋反應管<br>0.2 ml PCR Tubes with Flat Caps | 版本<br>Version | 第二版<br>Version II        |

| 產品說明 Specification |   |                 |                        | 產品圖 Drawing |
|--------------------|---|-----------------|------------------------|-------------|
| 產品尺寸<br>Dimension  | 管口直徑 / 管身高度<br>Ø 5.28 ±0.01mm / Height 21.33 ±0.1mm   |                 |                        |             |
| 顏色<br>Color        | 透明<br>Clear   |                 |                        |             |
| 包裝<br>Packaging    | 500 個/袋, 6 袋/盒, 10 盒/箱<br>500 pcs/ bag, 6 bags/ box, 10 boxes/ carton   |                 |                        |             |
| 產品容量<br>Volume     | 0.23 ml / tube  | 儲存條件<br>Storage | 室溫<br>Room temperature |             |
| 原料<br>Material     | <ul style="list-style-type: none"><li>● 本產品原料符合美國 FDA、美國藥典第六級認證之聚丙烯</li><li>● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)</li><li>● 不添加殺菌劑(二 (2-羥乙基), 二甲基十二烷基銨)</li><li>● 不添加增塑劑(鄰苯二甲酸鹽)</li><li>● MB-P02-A is made of FDA and USP VI standard virgin Polypropylene</li><li>● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)</li><li>● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMA))</li><li>● NOT INCLUDING: Plasticizers (phthalates)</li></ul>   |                 |                        |             |
| 產品特色<br>Features   | <ul style="list-style-type: none"><li>● 霧面蓋設計, 易於標記</li><li>● 薄壁設計, 導熱效果佳</li><li>● 無 DNA 分解酶, RNA 分解酶, 人類 DNA, PCR 抑制劑, ATP, 內毒素</li><li>● 可高溫高壓滅菌</li><li>● 獨特防漏設計, 確保樣品蒸散率&lt;5%</li><li>● 適用大部分主流 PCR 機</li><li>● 小包裝設計, 避免交叉污染</li><li>● Easy labelling on frosted flat cap</li><li>● Thin walled tube for best heat transfer</li><li>● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin</li><li>● Autoclavable</li><li>● Evaporation rate below 5%</li><li>● Small package can avoid contamination by different researchers</li></ul> |                 |                        |             |

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段192號9樓之1  
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,  
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332  
傳真(Fax): +886-2-2677-2339  
官網(Site): [www.gunster.com.tw](http://www.gunster.com.tw)  
電郵(Email): [gunster@gunster.com.tw](mailto:gunster@gunster.com.tw)

*Eiffie Wu*

Product Manager

吳佩儒  
產品經理