

技術資料表 Technical Data Sheet



GUNSTER BIOTECH

~Best Choice for Life Sciences~

貨號 Product NO.	MB-P08-AD	發布日 Date	2023/6/20 20-Jun-2023
品名 Description	0.2ml PCR 凸蓋獨立八連排 0.2ml PCR 8-Strip Tube with Attached Domed Caps	版本 Version	第二版 Version II

產品說明 Specification				產品圖 Drawing
產品尺寸 Dimension	管口直徑 / 管身高度 Ø 5.3 ±0.1mm / Height 23.3 ±0.1mm			
顏色 Color	透明 Clear			
包裝 Packaging	12 個/袋, 10 袋/盒, 20 盒/箱 12 pcs/ bag, 10 bags/ box, 20 boxes/ carton			
產品容量 Volume	0.22 ml / tube	儲存條件 Storage	室溫 Room temperature	
原料 Material	<ul style="list-style-type: none">● 本產品原料符合美國 FDA、美國藥典第六級認證之聚丙烯● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)● 不添加殺菌劑(二(2-羥乙基)，二甲基十二烷基銨)● 不添加增塑劑(鄰苯二甲酸鹽)● MB-P08-AD is made of FDA and USP VI standard virgin Polypropylene● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMA))● NOT INCLUDING: Plasticizers (phthalates)			
產品特色 Features	<ul style="list-style-type: none">● 適用凸蓋式熱循環儀器● 獨立八連排(管蓋一體)，減少汙染● 薄壁設計，導熱效果佳● 無 DNA 分解酶，RNA 分解酶，人類 DNA，PCR 抑制劑，ATP，內毒素● 可高溫高壓滅菌● 獨特防漏設計，確保樣品蒸散率<5%● 適用大部分主流 PCR 機● 小包裝設計，避免交叉汙染● Domed cap could fit hollow hot lid thermal cyclers● Attached caps help avoid cross-contamination between different tubes● Thin walled tube for best heat transfer● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin● Autoclavable● Evaporation rate below 5%● Small package can avoid contamination by different researchers			

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段192號9樓之1
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332
傳真(Fax): +886-2-2677-2339
官網(Site): www.gunster.com.tw
電郵(Email): gunster@gunster.com.tw

Eiffie Wu

Product Manager

吳佩儒
產品經理