


# 技術資料表 Technical Data Sheet

貨號 Product NO.	MB-ASM	發布日 Date	2026/3/16 16-March-2026
品名 Description	鋁箔抗凍封膜 Aluminum Membrane	版本 Version	第六版 Version VI

產品說明 Specification			產品圖 Drawing
產品尺寸 Dimension	長度 142 ±0.5 毫米 / 寬度 79.3 ±0.5 毫米 / 厚度(不含離層紙) 0.06 ±0.01 毫米 Length 142 ±0.5 mm / Width 79.3 ±0.5 mm / Thickness (Backing+Adhesive layer ) 0.06 ±0.01 mm		
	顏色 Color	銀色 Silver	
包裝 Packaging	10 片/袋, 10 袋/盒, 20 盒/箱 10 pcs/ bag, 10 bags/ box, 20 boxes/ carton	儲存條件 Storage	室溫 Room temperature
原料 Material	<ul style="list-style-type: none"><li>● 背膠：鋁箔</li><li>● 膠面：壓克力膠</li><li>● 離型層：木纖維、聚乙烯、矽</li><li>● Backing Material: Aluminum Foil</li><li>● Adhesive Material: Acrylic Adhesive</li><li>● Liner Material: Wood Fiber, Polyethylene, Silicone</li></ul>		
產品特色 Features	<ul style="list-style-type: none"><li>● 適用溫度：-80℃ ~ 120℃</li><li>● 適用於 PCR 或貯存樣品</li><li>● 黏性極佳，確保樣品不揮發</li><li>● 好穿刺設計，方便樣品取用</li><li>● 無 DNA 分解酶，RNA 分解酶，人類 DNA，PCR 抑制劑，ATP，內毒素</li><li>● 小包裝設計，避免交叉污染</li><li>● 鋼板黏著力 13.7N/25mm，聚丙烯 11.7N/25mm</li><li>● Operating temperature -80℃ ~ 120℃</li><li>● For PCR or Storage use</li><li>● Exceptional adhesion ensures sample integrity and prevents evaporation</li><li>● Ergonomically designed for easy penetration, facilitating convenient sample retrieval</li><li>● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin</li><li>● Small packages can avoid contamination by different researchers</li><li>● Adhesion to steel 13.7N/25mm, adhesion to polypropylene 11.7N/25mm</li></ul>		

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段192號9樓之1  
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,  
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332  
傳真(Fax): +886-2-2677-2339  
官網(Site): [www.gunster.com.tw](http://www.gunster.com.tw)  
電郵(Email): [gunster@gunster.com.tw](mailto:gunster@gunster.com.tw)

*Eiffie Wu*

Product Manager

吳佩儒  
產品經理