


貨號 Product NO.	MB-QR4	發布日 Date	2023/6/20 20-Jun-2023
品名 Description	qPCR 四連排 qPCR 4-Strip Tubes and Caps For Rotor-Gene® Q	版本 Version	第二版 Version II

產品說明 Specification				產品圖 Drawing
產品尺寸 Dimension	管口直徑 / 管身高度 Ø 3.36 ±0.1mm / Height 27.75 ±0.1mm			
顏色 Color	透明管+霧面蓋 Clear tube with Frosted Cap			
包裝 Packaging	252 個/袋, 2 袋/盒, 20 盒/箱 252 pcs/ bag, 2 bags/ box, 20 boxes/ carton			
產品容量 Volume	0.16 ml / tube	儲存條件 Storage	室溫 Room temperature	
原料 Material	<ul style="list-style-type: none"><li>● 本產品原料符合美國 FDA、美國藥典第六級認證之聚丙烯</li><li>● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)</li><li>● 不添加殺菌劑(二 (2-羥乙基), 二甲基十二烷基銨)</li><li>● 不添加增塑劑(鄰苯二甲酸鹽)</li><li>● MB-QR4 is made of FDA and USP VI standard virgin Polypropylene</li><li>● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)</li><li>● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMDA))</li><li>● NOT INCLUDING: Plasticizers (phthalates)</li></ul>			
產品特色 Features	<ul style="list-style-type: none"><li>● 超透明設計</li><li>● 薄壁設計, 導熱效果佳</li><li>● 無 DNA 分解酶, RNA 分解酶, 人類 DNA, PCR 抑制劑, ATP, 內毒素</li><li>● 獨特防漏設計, 確保樣品蒸散率&lt;5%</li><li>● 小包裝設計, 避免交叉污染</li><li>● 產品均含管身與上蓋</li><li>● Clear tubes</li><li>● Thin walled tube for best heat transfer</li><li>● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin</li><li>● Evaporation rate below 5%</li><li>● Small package can avoid contamination by different researchers</li><li>● This item contains tubes and caps</li></ul>			

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段 192 號 9 樓之 1  
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,  
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332  
傳真(Fax): +886-2-2677-2339  
官網(Site): [www.gunster.com.tw](http://www.gunster.com.tw)  
電郵(Email): [gunster@gunster.com.tw](mailto:gunster@gunster.com.tw)

*Eiffie Wu*

Product Manager

吳佩儒  
產品經理