

技術資料表 Technical Data Sheet



GUNSTER BIOTECH

~Best Choice for Life Sciences~

貨號 Product NO.	MB-384	發布日 Date	2023/6/20 20-Jun-2023
品名 Description	qPCR 384 孔盤(半裙邊) qPCR 384 Plate (ABI & Bio-Rad type)	版本 Version	第二版 Version II

產品說明 Specification				產品圖 Drawing
產品尺寸 Dimension	管口直徑 / 管身高度 Ø 3.2 ±0.1mm / Height 9.7 ±0.1mm (exclude the skirt)			
顏色 Color	透明霧 Frosted			
包裝 Packaging	10 個/袋, 5 袋/盒, 5 盒/箱 10 pcs/ bag, 5 bags/ box, 5 boxes/ carton			
產品容量 Volume	0.04 ml / tube	儲存條件 Storage	室溫 Room temperature	
原料 Material	<ul style="list-style-type: none">● 本產品原料符合美國 FDA、美國藥典第六級認證之聚丙烯● 不添加滑劑(油酰胺、芥酸酰胺、硬脂酰胺)● 不添加殺菌劑(二 (2-羥乙基), 二甲基十二烷基銨)● 不添加增塑劑(鄰苯二甲酸鹽)● MB-384 is made of FDA and USP VI standard virgin Polypropylene● NOT INCLUDING: Slip agents (including oleamide, erucamide, stearamide)● NOT INCLUDING: Biocides (di(2-hydroxyethyl), methyl dodecyl ammonium salts (DiHEMA))● NOT INCLUDING: Plasticizers (phthalates)			
產品特色 Features	<ul style="list-style-type: none">● 缺角位置：A24● 全襯裙● 薄壁設計，導熱效果佳● 無 DNA 分解酶，RNA 分解酶，人類 DNA，PCR 抑制劑，ATP，內毒素● 適用於 ABI 及 Bio-Rad 384 孔定量 PCR 儀● 小包裝設計，避免交叉污染● Notch：A24 (one-notch)● Full Skirted● Thin walled tube for best heat transfer● Free of DNase, RNase, DNA, PCR Inhibitor, ATP, Endotoxin● Compatible with ABI and Bio-Rad 384 well thermal cyclers● Small package can avoid contamination by different researchers			

鋼士特生物科技有限公司 Gunster Biotech Co., Ltd.

221412 新北市汐止區大同路三段 192 號 9 樓之 1
9F-1, No.192, Sec. 3, Datong Rd., Xizhi Dist.,
New Taipei City 221412, Taiwan (R.O.C.)

電話(Tel): +886-2-2677-2332
傳真(Fax): +886-2-2677-2339
官網(Site): www.gunster.com.tw
電郵(Email): gunster@gunster.com.tw

Eiffie Wu

Product Manager

吳佩儒
產品經理